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Appendix ۲ – Summary Table of IFU updates

Category & Section	Old IFU/Addendum	Revision
<p>۴,۳. Using Endo Therapy accessories</p> <p>۴,۲ Using Endo Therapy accessories (all scopes in the above chart with an asterisk)</p> <p>High-frequency cauterization treatment: Warning</p>	<p>Do not perform high-frequency cauterization while supplying oxygen. This may result in combustion during cauterization.</p>	<p>- Consult the instructions as described in the instruction manuals provided by the high frequency cauterization equipment manufacturer when performing high-frequency cauterization treatment.</p> <p>-Use caution when performing high-frequency cauterization while supplying oxygen as this may potentially result in combustion.</p> <p>-Avoid applying excessive energy to one spot as this may potentially result in combustion and/or charring of the tissue.</p> <p>- Use suction when there is smoke to secure visibility on bronchoscope and remove potential combustible materials from the treatment area.</p> <p>- Based on recommendations of clinicians and published literature, the oxygen concentration typically needs to be lower than ۴۰٪ before performing high-frequency cauterization treatment. If oxygen concentration is too high, combustion may potentially occur which could cause patient burns.</p> <p>- Based on recommendation of clinicians and published literature, the energy output level typically needs to be lower than ۴۰ Watts when performing high-frequency cauterization treatment. If the output level is too high, the endoscope's and/or accessory's insulation may be damaged and could potentially cause operator and/or patient burns. Also, combustion may potentially occur which could cause patient burns.</p> <p>- Ensure that there is separation of >۴cm between the endoscope and the tracheal tube when performing the procedure. Otherwise, there is a risk of accidentally damaging the tracheal tube during the procedure, or of the fire generated by the airway ignition spreading to it, which could potentially cause harm to the patient or damage to the equipment.</p>
	<p>Set the electrosurgical unit to the minimum necessary output level. If the output</p>	<p>Removed and further clarified above</p>

	level is too high, the endoscope's and/or accessory's insulation may be damaged and cause operator and/or patient burns.	
<p>۴,۳. Using Endo Therapy accessories</p> <p>۴,۲ Using Endo Therapy accessories (all scopes in the above chart with an asterisk)</p> <p>Argon plasma coagulation (APC): Warning</p> <p>*This does not apply to BF-F۲۶.</p>	<p>-The argon gas itself is neither combustible nor a promoter of combustible substances, but the argon plasma is very hot and could ignite combustible substances. Flammable substances burn easily when argon is irradiated in the presence of combustible gas such as high-concentration or pure oxygen. Be sure to observe the following cautions.</p> <p>-Before and during APC, do not feed oxygen or other combustible gases and liquids into the tracheobronchial system.</p>	<p>-Consult the instructions as described in the instruction manuals provided by the Argon Plasma Coagulation (APC) equipment manufacturer when performing APC treatment.</p> <p>-Use caution when performing APC treatment while supplying oxygen as this may potentially result in combustion.</p> <p>-Avoid applying excessive energy to one spot as this may potentially result in combustion and/or charring of the tissue.</p> <p>-Use suction when there is smoke to secure visibility on bronchoscope and remove potential combustible materials from the treatment area.</p> <p>-Based on recommendation of clinicians and published literature, the oxygen concentration typically needs to be lower than ۴۰٪ before performing APC treatment. If the oxygen concentration is too high, combustion may potentially occur which could cause patient burns.</p> <p>-Based on recommendation of clinicians and published literature, the energy output level typically needs to be lower than ۴۰ Watts when performing APC treatment. If the output level is too high, the endoscope's and/or accessory's insulation may be damaged and could potentially cause operator and/or patient burns. Also, combustion may potentially occur which could cause patient burns.</p> <p>-Ensure that there is separation of >۴cm between the endoscope and the tracheal tube when performing the procedure. Otherwise, there is a risk of accidentally damaging the tracheal tube during the procedure, or of the fire generated by the airway ignition spreading to it, which could potentially cause harm to the patient or damage to the equipment.</p>

<p>۴,۳. Using Endo Therapy accessories</p> <p>۴,۲ Using Endo Therapy accessories (all scopes in the above chart with an asterisk)</p> <p>Laser cauterization: Warning</p>	<p>-Do not perform laser cauterization while supplying oxygen. This may result in combustion during cauterization.</p> <p>-To avoid patient injury, burns, bleeding, perforation and/or damage to the endoscope, never emit laser radiation before confirming that an appropriate distance between the target and the endoscope's distal end is maintained and the tip of the laser probe is surely in the correct position in the endoscopic image.</p>	<p>-Olympus has only determined the compatibility of this endoscope with a Nd:YAG laser system. Therefore, it is only recommended to use this endoscope with a Nd:YAG laser system.</p> <p>-Consult the instructions as described in the instruction manuals provided by the laser equipment manufacturer when performing laser cauterization.</p> <p>-Use caution when performing laser cauterization while supplying oxygen as this may potentially result in combustion.</p> <p>-Avoid applying excessive energy to one spot as this may potentially result in combustion and/or charring of the tissue.</p> <p>-Use suction when there is smoke to secure visibility on bronchoscope and remove potential combustible materials from the treatment area.</p> <p>-Based on recommendation of clinicians and published literature, the oxygen concentration typically needs to be lower than ۴۰٪ before performing laser cauterization. If the oxygen concentration is too high, combustion may potentially occur which could cause patient burns.</p> <p>-Based on recommendation of clinicians and published literature, the energy output level typically needs to be lower than ۴۰ Watts when performing laser cauterization. If the output level is too high, combustion may potentially occur which could cause patient burns.</p> <p>-Ensure that there is separation of >۴cm between the endoscope and the tracheal tube when performing the procedure. Otherwise, there is a risk of accidentally damaging the tracheal tube during the procedure, or of the fire generated by the airway ignition spreading to it, which could potentially cause harm to the patient or damage to the equipment.</p> <p>-To avoid patient injury, burns, bleeding, perforation and/or damage to the endoscope, never emit laser radiation before confirming that the</p>
---	--	---

		endoscope's distal end is away from the target and the tip of the laser probe is surely in the correct position in the endoscopic image.
۲,۲ Specifications ۲,۳ Specifications (BF-F۲۶۰, BF-PE۲ & BF-TE۲)	۲.۲۳ Addendum: Only Nd:YAG laser or ۸۱۰ nm diode lasers may be used with Olympus laser compatible bronchoscopes. Olympus has not evaluated any other lasers for compatibility with the indicated bronchoscope models.	Compatible; Nd:YAG laser system only

High-frequency cauterization treatment

Endoscope model: all models except BF-XP190

WARNING

- Never use the high-frequency EndoTherapy accessories for the BF-XP190 because the distal end of the BF-XP190 is not insulated. Patient injury may result.
- Consult the instructions as described in the instruction manuals provided by the high-frequency cauterization equipment manufacturer when performing high-frequency cauterization treatment.
- Use caution when performing high-frequency cauterization while supplying oxygen as this may potentially result in combustion.
- Avoid applying excessive energy to one spot as this may potentially result in combustion and/or charring of the tissue.
- Use suction when there is smoke to secure visibility on bronchoscope and remove potential combustible materials from the treatment area.
- Based on recommendation of clinicians and published literature, the oxygen concentration typically needs to be lower than 40% before performing high-frequency cauterization treatment. If the oxygen concentration is too high, combustion may potentially occur which could cause patient burns.
- Based on recommendation of clinicians and published literature, the energy output level typically needs to be lower than 40 Watts when performing high-frequency cauterization treatment. If the output level is too high, the endoscope's and/or accessory's insulation may be damaged and could potentially cause operator and/or patient burns. Also, combustion may potentially occur which could cause patient burns.
- Ensure that there is separation of >4cm between the endoscope and the tracheal tube when performing the procedure. Otherwise, there is a risk of accidentally damaging the tracheal tube during the procedure, or of the fire generated by the airway ignition spreading to it, which could potentially cause harm to the patient or damage to the equipment.
- Always confirm that the tissue is an appropriate distance away from the distal end of the endoscope. If high-frequency cauterization is performed when the distal end of the endoscope contacts the tissue, patient injury, burns, bleeding, perforation, and equipment damage may occur.

WARNING

- Ensure that there is separation of >4cm between the endoscope and the tracheal tube when performing the procedure. Otherwise, there is a risk of accidentally damaging the tracheal tube during the procedure, or of the fire generated by the airway ignition spreading to it, which could potentially cause harm to the patient or damage to the equipment.
- If it is required to activate APC for a few seconds, change between APC and oxygen feeding.
- Keep the distal end of the APC probe in the endoscopic image before and during activation. Never activate APC in a position you cannot observe.
- Do not perform APC without gloves. Operator injury can result.
- Always confirm that the tissue is an appropriate distance away from the distal end of the endoscope. If APC is performed with the distal end of the endoscope in contact with the tissue, patient injury may occur.
- Make sure that the distal end of the APC probe always lies more than 10 mm from the endoscope's distal end (see Figure 4.9). The protrusion by 10 mm or more can be identified when the first black ring on the APC probe's distal end is visible in the endoscopic image. Otherwise, the treated region cannot be irradiated correctly and the endoscope may be damaged. Using a damaged endoscope may cause patient injury.


Figure 4.9

- Make sure that the patient plate is properly attached to patient's skin. Improper connection between the patient plate and patient's skin surface may cause burns. For further details on patient plates, refer to the patient plate's instruction manual.
- Do not bring the distal end of the APC probe close to a metal stent. The tissue around the metal stent may be burned.
- Before performing APC, inspect the surface of the endoscope for any dents, bulges or other irregularities. The endoscope may be damaged and cause patient and/or operator burns.

Argon plasma coagulation (APC)

Endoscope model: all models except BF-XP190

WARNING

- Consult the instructions as described in the instruction manuals provided by the Argon Plasma Coagulation (APC) equipment manufacturer when performing APC treatment.
- Use caution when performing APC treatment while supplying oxygen as this may potentially result in combustion.
- Avoid applying excessive energy to one spot as this may potentially result in combustion and/or charring of the tissue.
- Use suction when there is smoke to secure visibility on bronchoscope and remove potential combustible materials from the treatment area.
- Based on recommendation of clinicians and published literature, the oxygen concentration typically needs to be lower than 40% before performing APC treatment. If the oxygen concentration is too high, combustion may potentially occur which could cause patient burns.
- Based on recommendation of clinicians and published literature, the energy output level typically needs to be lower than 40 Watts when performing APC treatment. If the output level is too high, the endoscope's and/or accessory's insulation may be damaged and could potentially cause operator and/or patient burns. Also, combustion may potentially occur which could cause patient burns.

WARNING

- Always confirm that the electrode section of the electrosurgical accessory is at an appropriate distance from the distal end of the endoscope. Confirm that the entire green marking (in case of WLI observation mode) at the distal tip of the electrosurgical accessory can be observed on the endoscopic image. If the electrode is used when it is too close to the distal end of the endoscope, the endoscope and/or ancillary equipment may be damaged. Patient injury, burns, bleeding, perforation, and/or equipment damage may result.


Figure 4.8

- Be sure to contact the electrode section of the high-frequency EndoTherapy accessory with tissue when performing high-frequency cauterization treatment to prevent equipment damage and operator and/or patient burns.
- Improper connection between the patient plate and patient's skin surface may cause burns. For further details on the patient plate, refer to the patient plate's instruction manual.
- Do not perform high-frequency cauterization without gloves. Operator injury can result.
- Before performing high-frequency cauterization, inspect the surface of the endoscope for any dents, bulges, or other irregularities. Otherwise, patient injury, burns, bleeding, perforation, and/or equipment damage may result.
- When performing high-frequency cauterization, do not use the electrosurgical unit's SPRAY coagulation mode. The endoscope may be damaged, and it can cause patient and/or operator burns.

■ Laser cauterization

Endoscope model: all models except BF-P190 and BF-XP190

WARNING

- The BF-P190 and the BF-XP190 are incompatible with laser cauterization. Performing laser cauterization may cause patient injury and/or equipment damage.
- Olympus has only determined the compatibility of this endoscope with a Nd YAG laser system. Therefore, it is only recommended to use this endoscope with a Nd YAG laser system.
- Consult the instructions as described in the instruction manuals provided by the laser equipment manufacturer when performing laser cauterization.
- Use caution when performing laser cauterization while supplying oxygen as this may potentially result in combustion.
- Avoid applying excessive energy to one spot as this may potentially result in combustion and/or charring of the tissue.
- Use suction when there is smoke to secure visibility on bronchoscope and remove potential combustible materials from the treatment area.
- Based on recommendation of clinicians and published literature, the oxygen concentration typically needs to be lower than 40% before performing laser cauterization. If the oxygen concentration is too high, combustion may potentially occur which could cause patient burns.

WARNING

- Based on recommendation of clinicians and published literature, the energy output level typically needs to be lower than 40 Watts when performing laser cauterization. If the output level is too high, combustion may potentially occur which could cause patient burns.
- Ensure that there is separation of >4cm between the endoscope and the tracheal tube when performing the procedure. Otherwise, there is a risk of accidentally damaging the tracheal tube during the procedure, or of the fire generated by the airway ignition spreading to it, which could potentially cause harm to the patient or damage to the equipment.
- To avoid patient injury, burns, bleeding, perforation and/or damage to the endoscope, never emit laser radiation before confirming that the endoscope's distal end is away from the target and the tip of the laser probe is surely in the correct position in the endoscopic image.
- Always wear protective eyewear when performing laser cauterization. Otherwise, operator injury may occur.

OLYMPUS URGENT FIELD SAFETY NOTICE
FY26-901 BF scopes Global FCA Announcement- Bronchoscopes

Name & Address of Hospital/Medical Facility:

Department:

Inventory information (Model & Serial Number(s) of Bronchoscopes)

Model	Serial Number

Date:

Dear Sir or Madam

I herewith confirm the receipt of your Field Safety Notice.

Further I confirm that have transferred the content of attached FSN to all affected departments on witch this action has an impact. I understand the necessity to follow the steps.

Name (Signature)

Position:

Please scan/ email your complete paper form