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Appendix B. Intrinsic Resistance

Intrinsic resistance is defined as inherent or innate (not acquired) antimicrobial resistance, which is reflected in wild-type antimicrobial patterns of all or almost all representatives of a species. Intrinsic resistance is so common that susceptibility testing is unnecessary. For example, *Citrobacter* species are intrinsically resistant to ampicillin.

These tables can be helpful in at least three ways: 1) they provide a way to evaluate the accuracy of testing methods; 2) they aid in the recognition of common phenotypes; and 3) they can assist with verification of cumulative antimicrobial susceptibility test data. In the tables, an "R" occurring with an organism-antimicrobial agent combination means that strains should test resistant. A small percentage (1% to 3%) may appear susceptible due to method variation, mutation, or low levels of resistance expression.

A "susceptible" result should be viewed with caution. Ensure antimicrobial susceptibility test results and identification are accurate and reproducible. See Appendix A, footnote "a."

B1. Enterobacteriaceae

Antimicrobial Agent Organism	Ampicillin	Amoxicillin-clavulanate	Ampicillin-sulbactam	Piperacillin	Ticarcillin	Cephalosporin I: Cefazolin, Cephalothin	Cephamycins: Cefoxitin, Cefotetan	Cephalosporin II: Cefuroxime	Imipenem	Tetracyclines	Tigecycline	Nitrofurantoin	Polymyxin B Colistin	Aminoglycosides
<i>Citrobacter freundii</i>	R	R	R			R	R	R						
<i>Citrobacter koseri</i>	R			R	R									
<i>Enterobacter aerogenes</i>	R	R	R			R	R	R						
<i>Enterobacter cloacae</i> complex	R	R	R			R	R	R						
<i>Escherichia coli</i>	There is no intrinsic resistance to β -lactams in this organism.													
<i>Escherichia hermannii</i>	R				R									
<i>Hafnia alvei</i>	R	R	R			R	R							
<i>Klebsiella pneumoniae</i>	R				R									
<i>Morganella morganii</i>	R	R				R		R	*		R	R	R	
<i>Proteus mirabilis</i>	There is no intrinsic resistance to penicillins and cephalosporins in this organism.									*	R	R	R	R
<i>Proteus penneri</i>	R					R		R	*	R	R	R	R	R
<i>Proteus vulgaris</i>	R					R		R	*	R	R	R	R	R
<i>Providencia rettgeri</i>	R	R				R			*	R	R	R	R	R
<i>Providencia stuartii</i>	R	R				R			*	R	R	R	R	†
<i>Salmonella</i> and <i>Shigella</i> spp.	There is no intrinsic resistance to β -lactams in these organisms; refer to WARNING below for reporting.													
<i>Serratia marcescens</i>	R	R	R			R	R	R					R	R
<i>Yersinia enterocolitica</i>	R	R			R	R								

Appendix B. (Continued)

B1. (Continued)

WARNING: For *Salmonella* spp. and *Shigella* spp., aminoglycosides, first- and second-generation cephalosporins, and cephamycins may appear active *in vitro*, but are not effective clinically and should not be reported as susceptible.

* *Proteus* species, *Providencia* species, and *Morganella* species may have elevated minimal inhibitory concentrations to imipenem by mechanisms other than by production of carbapenemases. Isolates that test as susceptible should be reported as susceptible.

† *Providencia stuartii* should be considered resistant to gentamicin, netilmicin, and tobramycin but not intrinsically resistant to amikacin.

NOTE 1: Cephalosporins III, cefepime, aztreonam, ticarcillin-clavulanate, piperacillin-tazobactam, and the carbapenems are not listed, because there is no intrinsic resistance in *Enterobacteriaceae*.

NOTE 2: *Enterobacteriaceae* are also intrinsically resistant to clindamycin, daptomycin, fusidic acid, glycopeptides (vancomycin, teicoplanin), lipoglycopeptides (oritavancin, telavancin), linezolid, tedizolid, quinupristin-dalfopristin, rifampin, and macrolides (erythromycin, clarithromycin, and azithromycin). However, there are some exceptions with macrolides (ie, *Salmonella* and *Shigella* spp. with azithromycin).

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Table 4C. Disk Diffusion: Reference Guide to Quality Control Frequency

15-Replicate (3 x 5 Day) Plan: Acceptance Criteria and Recommended Action*

Number Out of Range With Initial Testing (based on 15 replicates)	Conclusion From Initial Testing (based on 15 replicates)	Number Out of Range After Repeat Testing (based on all 30 replicates)	Conclusion After Repeat Testing
0-1	Plan is successful. Convert to weekly QC testing.	N/A	N/A
2-3	Test another 3 replicates for 5 days.	2-3	Plan is successful. Convert to weekly QC testing.
≥4	Plan fails. Investigate and take corrective action as appropriate. Continue QC each test day.	≥4	Plan fails. Investigate and take corrective action as appropriate. Continue QC each test day.

* Assess each QC strain/antimicrobial agent combination separately.

Abbreviations: N/A, not applicable; QC, quality control.

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
Table 4D. Disk Diffusion: Troubleshooting Guide

This table provides guidance for troubleshooting and corrective action for out-of-range QC, primarily using antimicrobial susceptibility tests with MHA. Refer to M02-A12 (disk diffusion), Chapter 4, Quality Control and Quality Assurance for additional information. Out-of-range QC tests should first be repeated. If the issue is unresolved, this troubleshooting guide should be consulted regarding additional suggestions for troubleshooting out-of-range QC results. In addition, if unresolved, manufacturers should be notified of potential product problems.

General Comment

(1) QC organism maintenance: Avoid repeated subcultures. Retrieve new QC strain from stock. If using lyophilized strains, follow the maintenance recommendations of the manufacturer. Store *E. coli* ATCC® 35218 and *K. pneumoniae* ATCC® 700603 stock cultures at -60°C or below and prepare working cultures weekly (refer to M02-A12, Subchapter 4.4).

Antimicrobial Agent	QC Strain	Observation	Probable Cause	Comments/Action
Aminoglycosides	Any	Zone too small	pH of media too low	Acceptable pH range = 7.2-7.4 Avoid CO ₂ incubation, which lowers pH.
Aminoglycosides	Any	Zone too large	pH of media too high	Acceptable pH range = 7.2-7.4
Aminoglycosides	<i>P. aeruginosa</i> ATCC® 27853	Zone too small	Ca ⁺⁺ and/or Mg ⁺⁺ content too high	Use alternative lot of media.
Aminoglycosides	<i>P. aeruginosa</i> ATCC® 27853	Zone too large	Ca ⁺⁺ and/or Mg ⁺⁺ content too low	Use alternative lot of media.
Amoxicillin-clavulanate	<i>E. coli</i> ATCC® 35218	Zone too small	Clavulanate is labile. Disk has lost potency.	Use alternative lot of disks. Check storage conditions and package integrity.
Ampicillin	<i>E. coli</i> ATCC® 35218	Zone too large (should be no zone—resistant)	Spontaneous loss of the plasmid encoding the β-lactamase	See general comment (1) on QC organism maintenance.
β-Lactam group	Any	Zone initially acceptable, but decreases and possibly out of range over time	Disk has lost potency.	Use alternative lot of disks. Check storage conditions and package integrity. Imipenem, clavulanate, and cefaclor are especially labile.
Aztreonam Cefotaxime Cefpodoxime Ceftazidime Ceftriaxone	<i>K. pneumoniae</i> ATCC® 700603	Zone too large	Spontaneous loss of the plasmid encoding the β-lactamase	See general comment (1) on QC organism maintenance.
Cefotaxime-clavulanate Ceftazidime-clavulanate	<i>K. pneumoniae</i> ATCC® 700603	Negative ESBL test	Spontaneous loss of the plasmid encoding the β-lactamase	See general comment (1) on QC organism maintenance.
Penicillins	Any	Zone too large	pH of media too low	Acceptable pH range = 7.2-7.4 Avoid CO ₂ incubation, which lowers pH.
Penicillins	Any	Zone too small	pH of media too high	Acceptable pH range = 7.2-7.4
Carbencillin	<i>P. aeruginosa</i> ATCC® 27853	Zone too small	QC strain develops resistance after repeated subculture.	See general comment (1) on QC organism maintenance.
Ticarcillin-clavulanate	<i>E. coli</i> ATCC® 35218	Zone too small	Clavulanate is labile. Disk has lost potency.	Use alternative lot of disks. Check storage conditions and package integrity.

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
Table 4D. (Continued)

Antimicrobial Agent	QC Strain	Observation	Probable Cause	Comments/Action
Clindamycin	<i>S. aureus</i> ATCC® 25923	Zone too small	pH of media too low	Acceptable pH range = 7.2–7.4 Avoid CO ₂ incubation, which lowers pH.
Clindamycin	<i>S. aureus</i> ATCC® 25923	Zone too large	pH of media too high	Acceptable pH range = 7.2–7.4
Macrolides	<i>S. aureus</i> ATCC® 25923	Zone too small	pH of media too low	Acceptable pH range = 7.2–7.4 Avoid CO ₂ incubation, which lowers pH.
Macrolides	<i>S. aureus</i> ATCC® 25923	Zone too large	pH of media too high	Acceptable pH range = 7.2–7.4
Quinolones	Any	Zone too small	pH of media too low	Acceptable pH range = 7.2–7.4 Avoid CO ₂ incubation, which lowers pH.
Quinolones	Any	Zone too large	pH of media too high	Acceptable pH range = 7.2–7.4
Tetracyclines	Any	Zone too large	pH of media too low	Acceptable pH range = 7.2–7.4 Avoid CO ₂ incubation, which lowers pH.
Tetracyclines	Any	Zone too small	pH of media too high	Acceptable pH range = 7.2–7.4
Tetracyclines	Any	Zone too small	Ca ⁺⁺ and/or Mg ⁺⁺ content too high	Use alternative lot of media.
Tetracyclines	Any	Zone too large	Ca ⁺⁺ and/or Mg ⁺⁺ content too low	Use alternative lot of media.
Sulfonamides Trimethoprim Trimethoprim-sulfamethoxazole	<i>E. faecalis</i> ATCC® 29212	Zone ≤ 20 mm	Media too high in thymidine content	Use alternative lot of media.
Various	Various	Zone too small	Contamination Use of magnification to read zones	Measure zone edge with visible growth detected with unaided eye. Subculture to determine purity and repeat if necessary.
Various	Any	Many zones too large	Inoculum too light Error in inoculum preparation Media depth too thin MHA nutritionally unacceptable	Repeat using McFarland 0.5 turbidity standard or standardizing device. Check expiration date and proper storage if using barium sulfate or latex standards. Use agar with depth approximately 4 mm. Recheck alternate lots of MHA.
Various	Any	Many zones too small	Inoculum too heavy Error in inoculum preparation Media depth too thick MHA nutritionally unacceptable	Repeat using McFarland 0.5 turbidity standard or standardizing device. Check expiration date and proper storage if using barium sulfate or latex standards. Use agar with depth approximately 4 mm. Recheck alternate lots of MHA.
Various	Any	One or more zones too small or too large	Measurement error Transcription error Random defective disk Disk not pressed firmly against agar	Recheck readings for measurement or transcription errors. Retest. If retest results are out of range and no errors are detected, initiate corrective action.
Various	Various	Zone too large	Did not include lighter growth in zone measurement (eg, double zone, fuzzy zone edge)	Measure zone edge with visible growth detected with unaided eye.

Table 4D. (Continued)

Antimicrobial Agent	QC Strain	Observation	Probable Cause	Comments/Action
Various	<i>S. pneumoniae</i> ATCC® 49619	Zones too large. Lawn of growth scanty.	Inoculum source plate too old and contains too many nonviable cells. Plate used to prepare inoculum should be 18–20 hours.	Subculture QC strain and repeat QC test or retrieve new QC strain from stock.
Various	Any	One QC strain is out of range, but other QC organism(s) are in range with the same antimicrobial agent.	One QC organism may be a better indicator of a QC problem.	Retest this strain to confirm reproducibility of acceptable results. Evaluate with alternative strains with known MICs. Initiate corrective action with problem QC strain/antimicrobial agents.
Various	Any	Two QC strains are out of range with the same antimicrobial agent.	Indicates a problem with the disk	Use alternative lot of disks. Check storage conditions and package integrity.
Various	Any	Zones overlap.	Too many disks per plate	Place no more than 12 disks on a 150-mm plate and 5 disks on a 100-mm plate; for some fastidious bacteria that produce large zones, use fewer.

Abbreviations: ATCC®, American Type Culture Collection; ESBL, extended-spectrum β-lactamase; MHA, Mueller-Hinton agar; MIC, minimal inhibitory concentration; QC, quality control.

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
Scenario #1:

Ampicillin *E. coli* ATCC® 25922; acceptable range: 15 to 22 mm

Week	Day	Lot Number (Disks)	Lot Number (MHA)	Result	Action
1	1	3564	16481	18	
2	1	3564	16481	19	
3	1	3564	16481	18	
4	1	3564	16481	19	
5	1	3564	16481	14	Out of range. Repeat QC same day.
5	2	3564	16481	17	In range. Five acceptable in-range QC tests for <i>E. coli</i> ATCC® 25922 with ampicillin disks lot 3564 and MHA lot 16481. Resume weekly QC testing.

Abbreviations: ATCC®, American Type Culture Collection; MHA, Mueller-Hinton agar; QC, quality control.

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Scenario #2:

Ampicillin *E. coli* ATCC[®] 25922; acceptable range: 15 to 22 mm

Week	Day	Lot Number (Disks)	Lot Number (MHA)	Result	Action
1	1	9661	16922	18	
2	1	9661	16922	19	
3	1	9661	16922	14	Out of range. Repeat QC same day.
3	2	9661	16922	18	In range. Three acceptable in-range QC tests for <i>E. coli</i> ATCC [®] 25922 with ampicillin disks lot 9661 and MHA lot 16922. Repeat QC 2 more consecutive days.
3	3	9661	16922	18	In range.
3	4	9661	16922	17	In range. Five acceptable in-range QC tests for <i>E. coli</i> ATCC [®] 25922 with ampicillin disks lot 9661 and MHA lot 16922. Resume weekly QC testing.

Abbreviations: ATCC[®], American Type Culture Collection; MHA, Mueller-Hinton agar; QC, quality control.

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